


ANGIOGRAPHY PACK



Contents:

- A. 01 table drape 150x150cm
- B. 01 angio drape 250x340cm
- C. 03 reinforced gown in SMS size XL
- D. 04 towel
- E. 01 instrument cover 80x100cm
- F. 02 instrument cover 100x100cm
- G. 02 bowl 1000cc
- H. 12 gauze 10x10cm 12ply
- I. 01 sponge with handle
- J. 01 manifold 3ways OFF HP RA
- K. 01 guidewire 0.035 180cm J3mm
- L. 01 scalpel no. 11
- M. 01 HP line RA/LLF 150CM PUR
- N. 01 syringe 3 rings 12cc L.L. RA

Use Destination:

- Sterile patient protection and surgical field isolation and accessories during surgery angiography procedures with femoral/radial access